

ATP measurements

- 1) Collect samples (25 mid-second instar larvae; 5 adult flies).
- 2) Rinse several times with 1 ml cold PBS to remove all traces of food that might be attached to the outside of the animal. Larvae can be washed in a 1.5 ml microfuge tube, but adults should be rinsed in a 9-well glass plate. Transfer adult flies to a 1.5 ml microfuge tube. Carefully remove all liquid. For larvae, centrifuge at 3,000 x g and remove all remnants of PBS.
- 3) Prior to homogenizing the samples, prepare the ATP reaction mix (Molecular Probes ATP kit; A22066) by mixing 3.56 ml ddH₂O, 200 µl 20x reaction buffer, 40 µl 0.1 M DTT, 200 µl 10 mM D-luciferin, and 1 µl firefly luciferase. The resulting mix is sufficient for 40 reactions, and should be kept on ice and protected from light.
- 4) Animals are rapidly homogenized in 100 µl of homogenization buffer [6 M guanidine HCL, 100 mM Tris (pH 7.8), 4 mM EDTA] with a pellet pestle (Kontes; 749521-1500) on ice. A motor can be used to facilitate homogenization (Kontes; 749540-0000).
- 5) Remove 10 µl of homogenized sample to measure protein content with a Bradford assay. Keep samples on ice and do not heat-treat. Protein samples can be frozen and stored at -80°C for later analysis.
- 6) The remaining samples are boiled for 5 minutes and centrifuged for 3 min at maximum speed in a refrigerated tabletop centrifuge that has been prechilled to 4°C.
- 7) Transfer 10 µl of the supernatant into a 1.5 ml microfuge tube and dilute 1:10 with 90 µl dilution buffer [25 mM Tris (pH 7.8), 100 µM EDTA], then transfer 10 µl of the diluted supernatant to a second 1.5 ml tube and dilute 1:75 by adding 740 µl of dilution buffer (final dilution of 1:750). The diluted homogenate is centrifuged at 20,000 x g, and 10 µl of the supernatant is transferred to individual wells of a white, opaque 96 well plate (Corning; 3362).
- 8) Prepare a series of ATP standards by diluting the 5 mM ATP stock solution provided with the assay kit with ddH₂O (0, 0.01, 0.05, 0.1, 0.5, 1 µM). Add 10 µl of each ATP standard solution to the first row of the plate to provide a standard curve. ATP standards can be stored at -20°C for several weeks.
- 9) Start the assay by adding 100 µl of the luciferase reaction mix with a multichannel pipette and immediately begin measuring luminescence with a plate reader. A minimum of three sequential measurements should be made for the entire plate, and the results should be averaged.
- 10) Determine the ATP concentration by comparing the luminescence measurements for each sample to the ATP standard curve.