

Total protein measurement

This protocol is designed to measure total protein from fly homogenates prepared from triglyceride or carbohydrate measurements. Samples should be kept cold prior to this analysis.

- 1) Prepare a series of BSA standard stock solutions. Most restriction enzymes from NEB are supplied with a 10 mg/ml BSA stock, which can be used for this purpose. Also, BSA concentration can be verified using a spectrophotometer ($A_{280}=0.667$ for 1 mg/ml BSA). Generate a 0.5 mg/ml standard by diluting 5 μ l of 10 mg/ml BSA with 95 μ l of PBS (100 μ l final volume). Conduct a series of 2-fold serial dilutions into PBS (e.g. 50 μ l of 0.5 mg/ml diluted in 50 μ l PBS for a 0.25 mg/ml standard) to generate a 0.0625, 0.125, 0.25 and 0.5 mg/ml protein standards. Aliquots of BSA stock solution can be stored at -20°C .
- 2) Dilute the 10 μ l of homogenized sample between 1:10 to 1:20 in cold PBS. Centrifuge the samples at maximum speed for 5 minutes in a tabletop centrifuge to remove any precipitate.
- 3) Pipette 10 μ l of a PBS blank and the BSA dilution series into the first row of a clear-bottom 96 well plate. Transfer 10 μ l of supernatant from each homogenized fly sample into individual wells of the subsequent rows.
- 4) Add 200 μ l of Bio-Rad Protein Assay (diluted 1:5 with water) to each sample using a multichannel pipette.
- 5) Incubate the plate at room temperature for 5 min. Centrifuge the plate in an appropriate swing-bucket rotor to remove any air bubbles present in the samples. Use a plate reader to measure absorbance at 595 nm.
- 6) Determine the protein concentration for each sample based on the BSA standard curve.